

Ion mobility as an added dimension for toxicology screening

Petur W. Dalsgaard, Christian B. Mollerup, and Kristian Linnet

Section of Forensic Chemistry, University of Copenhagen, Denmark

Several forensic laboratories have started to use UHPLC-HR-MS for drug screening because of its advantages compared to GC-MS, HPLC-DAD and immunoassays. But compared to GC-MS, these UHPLC-HR-MS methods lack the access to large libraries, inter laboratory reproducible spectra, and standardized gradients. Most laboratories using UHPLC-HR-MS build in-house libraries or get libraries with up to 2500 compounds from the manufacturer of the instruments. HR-MS instruments can in theory screen for an unlimited number of compounds, but in practice it gives a lot of false positive findings and is very time consuming regarding to data analysis.

We set out to develop a practical UPLC-IMS-QTOF-MS drug screening by processing HDMS^E acquired data from forensic blood samples against a >4000 compounds library (Targeted and Suspect Screening). Also we set out to find a practical way to identify unknown peaks (Non-Targeted Screening).